

PlatR is a **tablet application** that:

- > Enables users to import or prepare their own pipetting plans
- > Navigates their pipetting on microtiter plates by illuminating the wells according to the plan
- > Tracks history and creates reports



PLATR PACKAGES

Depending on your lab's needs, you may choose between **three packages**.

Basic package has all you need to start using PlatR in the lab. **Advanced** and **Premium** packages contain additional accessories and modules and enable barcode reader integration into the pipetting workflow, full sample traceability, and automatic reports with timestamps.

	Basic	Advanced	Premium
Samsung Tablet	✓	✓	✓
PlatR application functionalities:			
> Custom pipetting plans			
> Excel plan import	✓		
> Pipetting sequence setup		✓	
> Well illumination guidance			✓
> Plan execution logs			
Cherry Picking module functionalities:			
> Quick access to validated plans			
> Sample ID scan and match			
> Automatic progress on the plate	30 days trial	✓	✓
> Full sample traceability			
> Automatic report with timestamps			
Diagnostic module functionalities:			
> Quick access to validated plan			
> Sample ID scan to allocate it on the plan			
> Multiple samples per well (sample pooling)	30 days trial	30 days trial	✓
> Full sample traceability			
> Automatic report with timestamps			
PlatR accessories			
Adhesive plate holder	1	3	5
Stand for adjustable inclination	×	✓	✓
Screen protector	×	✓	✓
Touch pen	×	✓	✓
Holder for PCR plates	×	✓	✓
Holder for ELISA strips	×	×	✓
Bluetooth footswitch	×	×	✓
Barcode reader	×	×	✓
Customer support	✓	✓	✓
Technical support (5 hours up to 1 year)	×	×	✓